

O w n e r s m a n u a l

LivingMS  
LivingAtaxia  
LivingCP  
**LivingSpinal**  
LivingDiabetic  
LivingAmputee  
LivingALS



**IntimateRider™**

*Bringing Couples Together*

Distributed Exclusively by Living Spinal

Congratulations on your purchase of the IntimateRider™ and/or RiderMate™. HealthPostures has designed the IntimateRider™ with you and your partner in mind. What you are about to experience will enhance your love life.

### Care and Maintenance

- Periodically check the IntimateRider™ and RiderMate™ welds. Make sure all bolts and fasteners are tightened securely.
- The IntimateRider™ and RiderMate™ covers are machine washable in cold water. Line or machine dry on low.

### Safety

- Be aware of pressure points. Use briefly (a few minutes) and inspect pressure points for redness.
- Discontinue use if redness is apparent. Continued use could result in pressure sores.
- Be aware of all moving parts and pinch points when folding.
- Do not use the IntimateRider™ or RiderMate™ in the shower, tub, hot tub, or pool.
- Do not stand or kneel on the IntimateRider™ or RiderMate™.
- This product is not intended to be motorized.
- IntimateRider™ or RiderMate™ not recommend for persons over 275lbs.



## How to Use

1. Transferring from (from wheelchair to the IntimateRider™.)
  - The first time you transfer make sure someone is with you.
  - Brace yourself by putting your hand on the middle of the seat.
  - DO NOT USE THE BACK FOR ASSISTANCE.
2. Slide on RiderMate™ Cover, make sure seam is on the side or underneath bottom of unit.
3. Once IntimateRider™ is set up, minimal upper body movement will cause motion.
4. Review the positions page on back cover for fun options.

*Be Creative, Have Fun, and Enjoy the Ride!*

---

## Accessories

### RiderMate™

Elevate your next position with RiderMate™. Face to face, back to front or whatever you can think of, you can achieve with RiderMate™. A perfect match with the IntimateRider™ or it can be part of something new on its own. RiderMate™ is that special extra something to spice up any occasion.

- Non Slip Leg Tips
- Machine Washable, Comfortable Fabric
- Fold for easy storage

### IntimateRider™ Slip Covers

With many fun colors and patterns to choose from, these soft, easy to wash and changeable slip covers add that personal touch to your IntimateRider™. Why choose just one. Keep your partner guessing what mood you will be in each night, will it be leopard for the wild side or soft blue for your gentle side.


Color Options for the IntimateRider™ Slip Covers & Rider Cushions

### IntimateRider™ Cushion

The specially designed IntimateRider cushion fits perfectly on IntimateRider™ to provide optimal seating comfort. This cushion slips on and off with ease and is machine washable.

### IntimateRider™ & RiderMate™ Carrying Bags

Make your getaways special by taking your IntimateRider™ and RiderMate™ with you in their discreetly designed carrying bags.



Order these accessories on-line now  
at: [www.livingspinal.com](http://www.livingspinal.com)  
or call (619)810-0010



## Positions

So many positions, where do you start? IntimateRider's easy sway gives you so many new ways to connect with your partner. See [www.intimaterider.com](http://www.intimaterider.com) for more fun positions.



125 E. Main Street, Belle Plaine, MN 56011  
Fax 952.873.3741 • Ph 952.873.6136  
[www.intimaterider.com](http://www.intimaterider.com)

# Warranty Information

1. Manufacturer will provide the only warranty associated with the product and provide parts, service, replacement units, and costs of shipping for the full period of the product warranty. Warranty and return details are stipulated below for all IntimateRider products:
2. New, unaltered and unused products are able to be returned within 30 days of receipt by the customer. The customer is responsible for the cost of return shipping as well as a 15% restocking fee deducted from the refund.
  - a. Definition of new, unaltered and unused condition is: without showing signs of wear or damage in any way. This includes damage as a result of poor/faulty return packaging
3. For hygienic reasons, any products that have been used even lightly are ineligible for return
4. Products that are delivered in dysfunctional condition due to shipping damage or some other obvious defect or damage must be photo documented and sent to Distributor within two weeks of receipt. It is the responsibility of the customer to open and inspect their product for obvious defects and properly report them before use of the product. If the customer uses the product despite obvious damages or defects to the product, products are ineligible for return or warranty replacement.
5. A warranty guaranteeing all IntimateRider products to be free from defects in materials or workmanship exists for a period of 90 days from date of shipment. Upholstery materials are warranted to be free from material defects and work at time of delivery, but because they can be damaged by normal use, they are warranted for 60 days. Only IntimateRider Division approved components validate this warranty. This warranty does not include labor charges incurred in replacement parts installation. Freight charges to the factory are at the expense of consumer to seller. No credit will be given for any repairs to damaged or defective items without prior approval of Manufacturer or Distributor. Manufacturer reserves the right to require that any such items be returned to the factory for inspection. This warranty is the only warranty on IntimateRider products and all other warranties, including implied warranties of marketability and fitness for a particular purpose are hereby disclaimed. Damage resulting from misuse, negligence, accident or alterations are excluded and will not be paid. Distributor or Manufacturer shall not be liable for any consequential or incidental damages whatsoever.
6. The warranty does not cover the following: light shipping damages, customer dissatisfaction with the functionality, use or appearance of the product, damage from installation of the product, lost or stolen packages, or products misplaced or lost by the customer.
7. In the event of a product or part recall, Manufacturer will send necessary replacements directly to dealers referencing PO numbers and it is then the responsibility of the dealer to distribute parts and/or products to the customers.



support@livingspinal.com  
Ph 619-810-0010  
www.livingspinal.com